

UiO IN-BIOS5000/9000

Basic file formats

Arvind Sundaram

Oct 29, 2024

Norwegian Sequencing Centre
OUS, Ullevål, Oslo



File formats

- ❖ *FASTA* - sequence data
- ❖ *FASTQ* - sequence data including quality (**SRA**)
- ❖ GTF / GFF, GFF3 - Gene feature & transfer format
- ❖ *SAM* / **BAM** - aligned / mapped data (**CRAM**)
- ❖ *VCF* / **BCF** - variant calling

Text format

Binary format

FASTA

- ❖ Represents a sequence - DNA / RNA / protein
- ❖ Two lines:
 - ❖ header starts '>'
 - ❖ Second line contains the sequence

header

```
>NM_001199335.1 Danio rerio toll-like receptor 21 (tlr21), mRNA  
CTAGACGAAGCTCTCACAGAGTGAAAGGCTTTGTAATGTTGGACTCAGTTTATAGGAATT  
ACATTTAAAAGTGGAAGTGAAAAACATGGCACACTCTGCATGCCACAACTGATACTAAA  
GGCCACATTCATCTGTCTCATAAACTTGCCTGCAGCTACAGTTTCAGAAGTTGCATAGA  
GATCCCAGATTCTAATCATACAATCTTTACATGTGTTAAAAGTTATGAACGAGACATAAC  
TGCGATTGTGAGTGATGTATTTCCCACTGCATTAAATCTTACAATCTCTCACA
```

Genomes
Transcriptomes

FASTQ

- ❖ Represents a sequence along with quality information
- ❖ Four lines:
 - ❖ header starts '@'
 - ❖ Second line contains the sequence
 - ❖ Quality header starts with '+'
 - ❖ Fourth line contains the quality for each base in the second line

Illumina read Phred 33

| Instrument ID | Flowcell ID | tile# | Y axis | index used | | | | | |
|---|-------------|------------|--------|------------|-------|-------|---|-------|--------|
| run# | lane# | X axis | read# | | | | | | |
| @J00146 | :32 | :HJF3HBBXX | :8 | :1101 | :1387 | :1121 | 1 | :N:0: | GCGATA |
| NGTTCACCTTGGAGACCTGCTGCGGTTATGAGTACGACCTGGCGTGAAAACCTATTCCTTCCA | | | | | | | | | |
| + | | | | | | | | | |
| #AAA<<F<F<7A<-7<JFJJJJJJJFJF<FAJA<FJAAF<<<JFJJ7FFJJJJJJAJJJJJ< | | | | | | | | | |

FASTQ quality encoding

[illegible]

S - Sanger Phred+33, raw reads typically (0, 40)
X - Solexa Solexa+64, raw reads typically (-5, 40)
I - Illumina 1.3+ Phred+64, raw reads typically (0, 40)
J - Illumina 1.5+ Phred+64, raw reads typically (3, 41)
with 0=unused, 1=unused, 2=Read Segment Quality Control Indicator (bold)
(Note: See discussion above).
L - Illumina 1.8+ Phred+33, raw reads typically (0, 41)

GFF/GTF

- ❖ GFF: General Feature Format
 - ❖ Current format is based on version 3 and called GFF3
 - ❖ format consists of one line per feature
 - ❖ each containing 9 columns of data, plus optional track definition lines.
- ❖ GTF: General Transfer Format is identical to GFF version 2.5.

GFF3

| seqid | source | chromosome coordinate system type | start | end | score | strand | phase | attributes | ID Name Parent Note *case-sensitive |
|-----------------|---------|---|-------|------|-------|--------|-------|---|---|
| ##gff-version 3 | | | | | | | | | |
| ctg123 | example | gene | 1050 | 9000 | . | + | . | ID=EDEN;Name=EDEN;Note=protein kinase | |
| ctg123 | example | mRNA | 1050 | 9000 | . | + | . | ID=EDEN.1;Parent=EDEN;Name=EDEN.1;Index=1 | |
| ctg123 | example | five_prime_UTR | 1050 | 1200 | . | + | . | Parent=EDEN.1 | |
| ctg123 | example | CDS | 1201 | 1500 | . | + | 0 | Parent=EDEN.1 | |
| ctg123 | example | CDS | 3000 | 3902 | . | + | 0 | Parent=EDEN.1 | |
| ctg123 | example | CDS | 5000 | 5500 | . | + | 0 | Parent=EDEN.1 | |
| ctg123 | example | CDS | 7000 | 7608 | . | + | 0 | Parent=EDEN.1 | |
| ctg123 | example | three_prime_UTR | 7609 | 9000 | . | + | . | Parent=EDEN.1 | |
| ctg123 | example | mRNA | 1050 | 9000 | . | + | . | ID=EDEN.2;Parent=EDEN;Name=EDEN.2;Index=1 | |
| ctg123 | example | five_prime_UTR | 1050 | 1200 | . | + | . | Parent=EDEN.2 | |
| ctg123 | example | CDS | 1201 | 1500 | . | + | 0 | Parent=EDEN.2 | |
| ctg123 | example | CDS | 5000 | 5500 | . | + | 0 | Parent=EDEN.2 | |
| ctg123 | example | CDS | 7000 | 7608 | . | + | 0 | Parent=EDEN.2 | |
| ctg123 | example | three_prime_UTR | 7609 | 9000 | . | + | . | Parent=EDEN.2 | |

BED format

- ❖ Browser Extensible Data
- ❖ 3-12 fields; but generally contains 6 or more

| chromosome # | chromStart | chromEnd | name | score | strand |
|--------------|------------|-----------|------|-------|--------|
| chr7 | 127471196 | 127472363 | Pos1 | 0 | + |
| Chr7 | 127472363 | 127473530 | Pos2 | 0 | + |
| chr7 | 127473530 | 127474697 | Pos3 | 0 | + |
| chr7 | 127474697 | 127475864 | Pos4 | 0 | + |
| chr7 | 127475864 | 127477031 | Neg1 | 0 | - |
| chr7 | 127477031 | 127478198 | Neg2 | 0 | - |

There are far more annotation formats but will stop here!!!

Chromosome coordinate systems

❖ 0-based

ACTGACTG

012345678

To represent the TGAC:

0-based inclusive: 2-5

1-based inclusive: 3-6

1-based exclusive: 3-7

Ensembl: 1-based inclusive

UCSC: 0-based start

1-based end

1-based for display

❖ 1-based

ACTGACTG

123456789

SAM, VCF, GFF, GTF: 1-based

BAM, BCF: 0-based

BED: 0-based exclusive

Most tools are aware of these differences

SAM/BAM

- ❖ SAM - Sequence Alignment/Map format
- ❖ BAM - same as/similar to SAM but in binary format
 - ❖ Most softwares and tools prefer BAM
- ❖ Header
 - (access using samtools view -H)
- ❖ Alignment records
 - (access using samtools view; use -h to include Header)

SAM/BAM: header

```
@HD VN:1.0GO:none SO:coordinate ← sort order
@SQ SN:chrM LN:16571
@SQ SN:chr1 LN:249250621
...
@SQ SN:chrX LN:155270560
@SQ SN:chrY LN:59373566
...
...
@RG ID:86-191PL:ILLUMINA LB:IL500 SM:86-191-1 ← Read group with platform,
@RG ID:Bsk010PL:ILLUMINA LB:IL501 SM:Bsk010-1 library and sample information
...
@RG ID:SDH023PL:ILLUMINA LB:IL508 SM:SDH023 ↙ Program analysis history
@PG ID:GATK IndelRealigner VN:2.0-39-gd091f72CL:knownAlleles=[ ]
targetIntervals=tmp.intervals.list LODThresholdForCleaning=5.0
consensusDeterminationModel=USE_READS entropyThreshold=0.15
maxReadsInMemory=150000 maxIsSizeForMovement=3000 maxPositionalMoveAllowed=200
maxConsensuses=30 maxReadsForConsensuses=120 maxReadsForRealignment=20000
noOriginalAlignmentTags=false nWayOut=null generate_nWayOut_md5s=false
check_early=false noPGTag=false keepPGTags=false indelsFileForDebugging=null
statisticsFileForDebugging=null SNPsFileForDebugging=null
@PG ID:bwaPN:bwaVN:0.6.2-r126
```


SAM/BAM: alignment records

| | | | | | | | | | | | | | | | | | |
|----|--|---|-----|---|------|---|-------|---|----|---|------|---|---|---|-------|---|------|
| 1 | HW-ST605:127:B0568ABXX:2:1201:10933:3739 | 2 | 147 | 3 | chr1 | 4 | 27675 | 5 | 60 | 6 | 101M | 7 | = | 8 | 27588 | 9 | -188 |
| 10 | TCATTTTATGGCCCCTTCTTCCTATATCTGGTAGCTTTTAAATGATGACCATGTAGATAATCTTTATTGTCCCTCTTTCAGCAGAC | | | | | | | | | | | | | | | | |
| 11 | =7;::;<=?<=BCCEFFEJFCEGGEFFDF?E@E@ADCACB>CCDCBACDCDDAB@@BCADDCBC@BCBB8@ABCCDCBDA@>:/ | | | | | | | | | | | | | | | | |
| 12 | RG:Z:86-191 | | | | | | | | | | | | | | | | |

| Col | Field | Type | Regexp/Range | Brief Description |
|-----|-------|--------|--|---------------------------------------|
| 1 | QNAME | String | [!-?A-~]{1,254} | Query template NAME |
| 2 | FLAG | Int | [0,2 ¹⁶ -1] | bitwise FLAG |
| 3 | RNAME | String | * [!-()+-<>-~][!-~]* | Reference sequence NAME |
| 4 | POS | Int | [0,2 ³¹ -1] | 1-based leftmost mapping POSition |
| 5 | MAPQ | Int | [0,2 ⁸ -1] | MAPping Quality |
| 6 | CIGAR | String | * ([0-9]+ | CIGAR string |
| 7 | RNEXT | String | * = [!-()+-<>-~][!-~]* | Ref. name of the mate/next read |
| 8 | PNEXT | Int | [0,2 ³¹ -1] | Position of the mate/next read |
| 9 | TLEN | Int | [-2 ³¹ +1,2 ³¹ -1] | observed Template LENght |
| 10 | SEQ | String | * [A-Za-z=.]+ | sequence SEQUENCE |
| 11 | QUAL | String | [!-~]+ | ACSII of Phred-scaled base QUALity+33 |

SAM/BAM

- ❖ To check FLAG info: <https://broadinstitute.github.io/picard/explain-flags.html>
- ❖ BAM/SAM files are produced by alignment/mapping softwares and other tools.
- ❖ QC of BAM/SAM files - Picard, Qualimap
- ❖ Manipulating and extracting information from varied file formats - SAMtools, BEDtools, BCFtools

VCF/BCF

- ❖ Variant call format
- ❖ BCF - same as/similar to VCF but in binary format
 - ❖ Most softwares and tools prefer BCF or VCF in zip format
- ❖ header and records

header

```
##fileformat=VCFv4.2
##source=myImputationProgramV3.1
##reference=file:///seq/references/1000GenomesPilot-NCBI36.fasta
##contig=<ID=20,length=62435964,assembly=B36,md5=xxxx,species="Homo sapiens",taxonomy=x>
...
##INFO=<ID=NS,Number=1,Type=Integer,Description="Number of Samples With Data">
##INFO=<ID=DP,Number=1,Type=Integer,Description="Total Depth">
...
##FILTER=<ID=q10,Description="Quality below 10">
...
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
...
#CHROM POS ID REF ALT QUAL FILTER INFO FORMAT NA00001 NA00002 NA00003
```

VCF/BCF: records

```
#CHROM POS ID REF ALT QUAL FILTER INFO FORMAT NA00001
20 14370 rs6054257 G A 29 PASS NS=3;DP=14;AF=0.5;DB;H2 GT:GQ:DP:HQ 0|0:48:1:51,51

20 17330 . T A 3 q10 NS=3;DP=11;AF=0.017 GT:GQ:DP:HQ 0|0:49:3:58,50

20 1110696 rs6040355 A G,T 67 PASS NS=2;DP=10;AF=0.333,0.667;AA=T;DB GT:GQ:DP:HQ 1|2:21:6:23,27
```

| | #CHRO | POS | ID | REF | ALT | QUAL | FILTER | INFO | FORMAT | SAMPLEs |
|-----------|-------|-----|----|-----|----------|------|--------|-------|--------|---------|
| SNP | 20 | 3 | . | C | G | . | PASS | DP=10 | | |
| Deletion | 20 | 2 | . | TC | C | . | PASS | DP=10 | | |
| Insertion | 20 | 2 | . | TC | TCA | . | PASS | DP=10 | | |
| Alleles | 20 | 2 | . | TC | TG,T | . | PASS | DP=10 | | |
| Alleles | 20 | 2 | . | TCG | TG,T,TCA | . | PASS | DP=10 | | |

+ Structural variants

Qualimap

- ❖ GUI, command line based
 - ❖ examines sequencing alignment data in SAM/BAM format
- ❖ requires Java, R
- ❖ Uses:
 - ❖ FastQC, Picard, SAMTools, NOISeq, Repitools and JProfiler
 - ❖ BAM QC
 - ❖ RNA-seq QC
 - ❖ Counts QC
 - ❖ Multi-sample BAM QC

SAMtools/BCFtools

❖ SAMtools

- ❖ manipulating SAM/BAM

❖ BCFtools

- ❖ manipulating VCF/BCF

❖ Part of HTSlib tool kit

```
-- indexing
  faidx    index/extract FASTA
  index    index alignment
-- editing
  calmd    recalculate MD/NM tags and '=' bases
  fixmate  fix mate information
  reheader replace BAM header
  rmdup    remove PCR duplicates
  targetcut cut fosmid regions (for fosmid pool only)
-- file operations
  bamshuf  shuffle and group alignments by name
  cat      concatenate BAMs
  merge    merge sorted alignments
  mpileup  multi-way pileup
  sort     sort alignment file
  split    splits a file by read group
  bam2fq   converts a BAM to a FASTQ
-- stats
  bedcov   read depth per BED region
  depth    compute the depth
  flagstat simple stats
  idxstats BAM index stats
  phase    phase heterozygotes
  stats    generate stats (former bamcheck)
-- viewing
  flags    explain BAM flags
  tview    text alignment viewer
  view     SAM<->BAM<->CRAM conversion
```

BEDtools

- ❖ Awesome tool to compare dataset in BED / GFF / VCF and BAM formats
- ❖ Very well documented and Highly recommended!!
- ❖ Example:

<http://bedtools.readthedocs.io/en/latest/content/tools/intersect.html>



IGV - Integrative Genomics Viewer




- ❖ View BAM, VCF, GTF files and more

- ❖ Ensembl browser

<http://ensembl.org/>

- ❖ UCSC browser

<http://genome.ucsc.edu/>



<http://software.broadinstitute.org/software/igv/>

Further reading

- ❖ GFF3, GTF - <http://gmod.org/wiki/GFF3>
- ❖ SAM/BAM format - <http://samtools.github.io/hts-specs/SAMv1.pdf>
- ❖ VCF/BCF format - <https://samtools.github.io/hts-specs/VCFv4.2.pdf>
- <http://vcftools.sourceforge.net/VCF-poster.pdf>
- ❖ IGV - <http://software.broadinstitute.org/software/igv/>

Databases

- ❖ Ensembl - <http://ensembl.org/>
- ❖ UCSC - <http://genome.ucsc.edu/>
- ❖ NCBI - <http://ncbi.nlm.nih.gov/>
 - ❖ RefSeq - <https://www.ncbi.nlm.nih.gov/refseq/>

There are too many databases containing specific information!!!