

UiO IN-BIOS5000/9000

Illumina Technology

Arvind Sundaram

Oct 28, 2024

Norwegian Sequencing Centre
OUS, Ullevål, Oslo



High Throughput Sequencing

- ❖ First generation - **past, present**
 - ❖ Up to 1 kb; high quality data; multiplexed
 - ❖ SANGER; Highly automated (ABI Sanger 3730xl)
- ❖ Second generation - **present**
 - ❖ Shorter reads; Massive parallelisation and real high throughput
 - ❖ Illumina, BGISEq, Ion-torrent, [454, Solid]
 - ❖ RNA is reverse-transcribed to cDNA before sequencing
- ❖ Third generation - **[present] future**
 - ❖ Long-read sequencing; Single-molecule sequencing (without amplification)
 - ❖ PacBio, Oxford Nanopore, [more in development]
 - ❖ Potential to sequence RNA directly


Illumina sequencers

- ❖ Second generation sequencing technique
- ❖ Sequencing-by-synthesis aka SBS
 - ❖ <https://www.youtube.com/watch?v=fCd6B5HRaZ8>
- ❖ Mass parallelisation and real high throughput

Library prep and sequencing


Fragment (DNA) sequenced: up to 800 bp

A. Library Preparation




NGS library is prepared by fragmenting a gDNA sample and ligating specialized adapters to both fragment ends.

B. Cluster Amplification




Library is loaded into a flow cell and the fragments are hybridized to the flow cell surface. Each bound fragment is amplified into a clonal cluster through bridge amplification.

C. Sequencing



Sequencing reagents, including fluorescently labeled nucleotides, are added and the first base is incorporated. The flow cell is imaged and the emission from each cluster is recorded. The emission wavelength and intensity are used to identify the base. This cycle is repeated "n" times to create a read length of "n" bases.



Sequencing by synthesis (SBS)






Old sequencers use 4 colors
Newer machines use 2 colors
iSeq uses 1 color

XLEAP-SBS (2022)

- Increased speed
- Greater fidelity
- Greater accuracy
- More robustness




Sequencing



Library prep and sequencing

Fragment (DNA) sequenced: up to 800 bp



Read type: Mate pair



Add adapters during library preparation

Multiplexing: single / dual index




Multiplexing: pooling



Read type: Single end



Read type: Paired end




Single/Paired-end/Mate pair sequencing


- ❖ Paired-end reads or mate pair reads are pairs of reads known to come from two close regions in the genome.
- ❖ They are located with an approximate fixed distance from each other.
- ❖ Typically paired ends are a ~100-500bp apart, while mate pairs are ~2-10kb apart
- ❖ Allows short reads to have a larger "effective" size
- ❖ Performed by sequencing fragments from both ends
- ❖ Often used with Illumina reads
 - ❖ Typically 2 x 150 bp separated by 300bp
- ❖ May also overlap (e.g. 2x250bp from 400bp fragments)

Paired-end sequencing

- ❖ Both ends of fragments will be sequenced
- ❖ Gives information on genomic distance between pairs of reads
- ❖ May be used to overcome some problems with short reads



Adapters and multiplexing



DRAGEN on-board/external



BCL Convert
DRAGEN ORA Compression
DRAGEN FASTQC + MultiQC
Whole Genome
Enrichment (Including Exome)
DRAGEN Amplicon
RNA
Single Cell RNA
NanoString GeoMx NGS
Methylation
Metagenomics
RNA Pathogen Detection
COVID
TSO 500 Portfolio
Imputation
ScATAC-Seq
PGx Star Allele Caller
Illumina Complete Long Reads
DRAGEN secondary analysis for RIPIP and UPIP



Cloud based analysis

Stand-alone server

Integrated on-board

- NextSeq 1000/2000
- NovaSeq X
- MiSeq i100

FPGA based-architecture

30 minute* end-to-end analysis

Now open-source (2024)

Illumina sequencers - benchtop



iSeq 100 System



MiniSeq System



MiSeq System



MiSeq i100 and MiSeq i100 Plus Systems^a

2024


XLEAP-SBS
DRAGEN

Flow cell	i1	Mid-output	Rapid	High-output	Nano	Micro	v2	v3	5M	25M	50M	100M
Output range	144 Mb–1.2 Gb	2.1–2.4 Gb	2 Gb	1.65–7.5 Gb	300–500 Mb	1.2 Gb	750 Mb–8.5 Gb	3.8–15 Gb	1.5–3 Gb	2.5–15 Gb	5–30 Gb	10–30 Gb
Single-end reads per run	4M	8M	20M	25M	1M	4M	15M	25M	5M	25M	50M	100M
Run time (hr) ^b	9–19	17	< 5	7–24	17–28	19	5.5–39	21–56	7–15	4–15	4–15	5–8
Maximum read length (bp)	2 × 150	2 × 150	1 × 100	2 × 150	2 × 250	2 × 150	2 × 250	2 × 300	2 × 300	2 × 300	2 × 300	2 × 150
Included data analysis	Local Run Manager	Local Run Manager			Local Run Manager				DRAGEN™ software			

a. The MiSeq i100 System supports the 5M and 25M flow cells only; the MiSeq i100 Plus System supports all four flow cells.

b. Listed run times are estimates.

Data output



Illumina sequencers - slightly bigger



NextSeq 550 System^a



NextSeq 1000 and NextSeq 2000 Systems


2020

XLEAP-SBS

DRAGEN

Flow cell	Mid-output	High-output	P1 ^b	P2 ^b	P3 ^c	P4 ^c
Output range	16–39 Gb	25–120 Gb	10–60 Gb	40–240 Gb	120–360 Gb	90–540 Gb
Single-end reads per run	130M	400M	100M	400M	1.2B	1.8B
Run time (hr)	15–26	11–29	8–34	12–42	18–40	12–44
Maximum read length (bp)	2 × 150	2 × 150	2 × 300	2 × 300	2 × 150	2 × 150
Included data analysis	Local Run Manager		Onboard DRAGEN secondary analysis			

Data output



Illumina sequencers - production scale



NovaSeq 6000 System



NovaSeq X System




NovaSeq X Plus System

2022
XLEAP-SBS
DRAGEN

2022
XLEAP-SBS
DRAGEN


Flow cell	SP	S1	S2	S4	1.5B	10B	25B	1.5B	10B	25B
Flow cells processed per run	1 or 2	1 or 2	1 or 2	1 or 2	1	1	1	1 or 2	1 or 2	1 or 2
Output range	65–800 Gb	134 Gb–1 Tb	333 Gb–2.5 Tb	280 Gb–6 Tb	165–500 Gb	1–3 Tb	8 Tb	165 Gb–1 Tb	1–6 Tb	8–16 Tb
Single-end reads per flow cell	800M	1.6B	4.1B	10B	1.6B	10B	26B	1.6B	10B	26B
Run time (hr)	13–38	13–25	16–36	< 44	17–23	18–25	~48	17–23	18–25	~48
Maximum read length (bp)	2 × 250	2 × 150	2 × 150	2 × 150	2 × 150	2 × 150	2 × 150	2 × 150	2 × 150	2 × 150
Included data analysis	—				Onboard DRAGEN secondary analysis					

Data output



Data output

Illumina sequencers



What can you sequence using Illumina

- ❖ DNA studies
 - ❖ Whole genome sequencing - short reads are a pitfall
 - ❖ Genome re-sequencing
 - ❖ Exomes and target re-sequencing...
 - ❖ ChiP seq and more...
- ❖ RNA studies
- ❖ modification studies
 - ❖ Methylation and more...



Applications

DRAGEN

DRAGEN

DRAGEN

iSeq™ 100 MiniSeq™ MiSeq™ MiSeq™ i100/i100 Plus NextSeq™ 550 NextSeq™ 1000/2000 NovaSeq™ 6000 NovaSeq™ X/X Plus

SMALL WHOLE-GENOME SEQUENCING

TARGETED GENE SEQUENCING

TARGETED GENE EXPRESSION PROFILING

16S METAGENOMIC SEQUENCING

EXOME SEQUENCING

TRANSCRIPTOME SEQUENCING

CELL-FREE SEQUENCING*

SINGLE-CELL OR SPATIAL PROFILING

SHOTGUN METAGENOMICS

METHYLATION ANALYSIS

CHROMATIN ANALYSIS*

LARGE WHOLE-GENOME SEQUENCING

- DNA
- RNA
- Epigenetics
- Other